



- The Lightest Bionic Hand at 700 grams
- Feel The difference between Soft and Hard Objects
- Learn to use Grippy in Under 30 Mins
- Pick up any shape of object with Adaptive Grips
- MMG Based Sensor for better Control than EMG
- IoT capability for seamless Evaluation and Analytics



Call us on:

+91-86683-72127



Muscle Pre-Diagnostic Tool for Objective Evaluation of Amputees





- Patient Information Management System
- Guided Process for Evaluation within 30 mins
- Seamless Measurements of Muscle Strength
- Identification of Peak Muscle Signals for Prosthesis
- Wireless Interface with Grippy for Rehabilitation



Call us on:

+91-86683-72127